

Product Datasheet

Human ATG7 ELISA Kit (orb405636)

Description	Human Ubiquitin-like modifier-activating enzyme ATG7(ATG7) ELISA kit
Reactivity	Human
Range	28 pg/mL-1800 pg/mL
Target	ATG7
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	O95352
Sensitivity	7 pg/mL
Expiration Date	Please enquire.