

Product Datasheet

Human ADAMTS9 ELISA Kit (orb405639)

Description	Human A disintegrin and metalloproteinase with thrombospondin motifs 9(ADAMTS9) ELISA kit
Reactivity	Human
Range	0.78 ng/mL-50 ng/mL
Target	ADAMTS9
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q9P2N4
Sensitivity	0.195 ng/mL
Expiration Date	Please enquire.