

Product Datasheet

Human EPHA2 ELISA Kit (orb405640)

Description	Human Ephrin type-A receptor 2 (EPHA2) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	EPHA2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P29317
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.