

Product Datasheet

Human C1QTNF3 ELISA Kit (orb405657)

Description	Human C1q/TNF-related protein-3, Cartonectin/CTRP3/CORS-26 ELISA Kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	C1QTNF3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates
Uniprot ID	Q9BXJ4
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.