

Product Datasheet

Human CCL7 ELISA Kit (orb405665)

Description	Human C-C motif chemokine 7 (CCL7/MCP3/SCYA6/SCYA7) ELISA kit
Reactivity	Human
Range	15.6 pg/mL-1000 pg/mL
Target	CCL7
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma,tissue homogenates
Uniprot ID	P80098
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.