

Product Datasheet

Human CD109 ELISA Kit (orb405670)

| | |
|--------------------------|--|
| Description | Human CD109 antigen(CD109) ELISA kit |
| Reactivity | Human |
| Range | 94 pg/mL-6000 pg/mL |
| Target | CD109 |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma, tissue homogenates, cell lysates |
| Uniprot ID | Q6YHK3 |
| Sensitivity | 23.4 pg/mL |
| Expiration Date | Please enquire. |