

Product Datasheet

Human CLEC7A ELISA Kit (orb405708)

Description	Human C-type lectin domain family 7 member A(CLEC7A) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	CLEC7A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9BXN2
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.