

Product Datasheet

Human DDIT4 ELISA Kit (orb405712)

Description	Human DNA damage-inducible transcript 4 protein(DDIT4) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	DDIT4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, tissue homogenates, cell lysates, urine, cerebrospinal fluid (CSF)
Uniprot ID	Q9NX09
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.