

Product Datasheet

Human TRIM63 ELISA Kit (orb405718)

Description	Human E3 ubiquitin-protein ligase TRIM63(TRIM63) ELISA kit
Reactivity	Human
Range	12.5 pg/mL-800 pg/mL
Target	TRIM63
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q969Q1
Sensitivity	3.1 pg/mL
Expiration Date	Please enquire.