

Product Datasheet

Human Gremlin 1 ELISA Kit (orb405739)

Description	Human Gremlin-1(GREM1) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	GREM1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
Uniprot ID	O60565
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.