

Product Datasheet

Human Collagen III ELISA Kit (orb405759)

Description	Human Collagen, type III, alpha 1 (COL3A1)ELISA Kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	COL3A1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P02461
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.