

Product Datasheet

Human SIRT3 ELISA Kit (orb405778)

Description	Human NAD-dependent deacetylase sirtuin-3, mitochondrial(SIRT3) ELISA kit
Reactivity	Human
Range	15.6 pg/mL-1000 pg/mL
Target	SIRT3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9NTG7
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.