

Product Datasheet

Human p53 ELISA Kit (orb405792)

Description	Human p53/tumor protein,p53/TP53 ELISA Kit
Reactivity	Human
Range	9.38 pg/mL-600 pg/mL
Target	TP53
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, cell culture supernates, urine, cerebrospinal fluid (CSF), tissue homogenates, cell lysates
Uniprot ID	P04637
Sensitivity	2.34 pg/mL
Expiration Date	Please enquire.