

## Product Datasheet

### Human P53 ELISA Kit (orb405793)

|                          |                              |
|--------------------------|------------------------------|
| <b>Description</b>       | Human P53 ELISA Kit          |
| <b>Reactivity</b>        | Human                        |
| <b>Range</b>             | 1.56 ng/mL-100 ng/mL         |
| <b>Target</b>            | P53 Ab                       |
| <b>Note</b>              | For research use only        |
| <b>Application notes</b> | Detection Wavelength: 450 nm |
| <b>Sample Types</b>      | serum, plasma                |
| <b>Sensitivity</b>       | 0.39 ng/mL                   |
| <b>Expiration Date</b>   | Please enquire.              |