

Product Datasheet

Human MDR1 ELISA Kit (orb405802)

Description	Human permeability glycoprotein,P-gp/ABCB1 ELISA Kit
Reactivity	Human
Range	3.12 ng/mL-200 ng/mL
Target	ABCB1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, ascitic fluid, tissue homogenates
Uniprot ID	P08183
Sensitivity	0.78 ng/mL
Expiration Date	Please enquire.