

Product Datasheet

Human ROCK2 ELISA Kit (orb405812)

Description	Human Rho-associated coiled-coil containing protein kinase 2(Rock-2) ELISA Kit
Reactivity	Human
Range	0.625 ng/mL-40 ng/mL
Target	ROCK2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates
Uniprot ID	075116
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.