

Product Datasheet

Human SMOC2 ELISA Kit (orb405841)

Description	Human SPARC-related modular calcium-binding protein 2(SMOC2) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	SMOC2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q9H3U7
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.