

Product Datasheet

Human TARDBP ELISA Kit (orb405847)

Description	Human TAR DNA-binding protein 43 (TARDBP/TDP43) ELISA kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Target	TARDBP
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, cerebrospinal fluid (CSF)
Uniprot ID	Q13148
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.