

Product Datasheet

Human KLF10 ELISA Kit (orb405850)

Description	Human TGF-beta-inducible early response gene-1,TIEG1 ELISA Kit
Reactivity	Human
Range	34.38 pg/mL-2200 pg/mL
Target	KLF10
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q13118
Sensitivity	8.6 pg/mL
Expiration Date	Please enquire.