

Product Datasheet

Human IMP3 ELISA Kit (orb405870)

Description	Human U3 small nucleolar ribonucleoprotein protein IMP3(IMP3) ELISA kit
Reactivity	Human
Range	39 pg/mL-2500 pg/mL
Target	IMP3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9NV31
Sensitivity	9.75 pg/mL
Expiration Date	Please enquire.