

Product Datasheet

Human alpha II-spectrin breakdown product SBDP145 ELISA Kit (orb405885)

Description	Human alpha II-spectrin breakdown product SBDP145 ELISA kit
Reactivity	Human
Range	0.625 ng/mL-40 ng/mL
Target	SBDP145
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cerebrospinal fluid (CSF)
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.