

## Product Datasheet

### Human ENO1 antibody (IgA) ELISA Kit (orb405900)

<b>Description</b>	Human alpha-enolase (ENO1) antibody (IgA) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	Request Information
<b>Target</b>	ENO1 Ab (IgA)
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma
<b>Sensitivity</b>	Request Information
<b>Expiration Date</b>	Please enquire.