

Product Datasheet

Human ENO1 ELISA Kit (orb405903)

Description	Human Alpha-enolase (ENO1/ENO1L1/MBPB1/MPB1) ELISA kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Target	ENO1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, tissue homogenates, plasma
Uniprot ID	P06733
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.