

Product Datasheet

Human ST6GAL1 ELISA Kit (orb405912)

Description	Human Beta-galactoside alpha-2,6-sialyltransferase 1(ST6GAL1) ELISA kit
Reactivity	Human
Range	3.9 ng/mL-250 ng/mL
Target	ST6GAL1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P15907
Sensitivity	0.98 ng/mL
Expiration Date	Please enquire.