

Product Datasheet

Human amyloid beta peptide 1-40 ELISA Kit (orb405915)

Description	Human amyloid beta peptide 1-40, A β 1-40 ELISA Kit
Reactivity	Human
Range	125 pg/mL-8000 pg/mL
Target	A β 1-40
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, tissue homogenates, cerebrospinal fluid (CSF)
Sensitivity	31.25 pg/mL
Expiration Date	Please enquire.