

Product Datasheet

Human amyloid beta peptide 1-42 ELISA Kit (orb405916)

Description	Human amyloid beta peptide 1-42, A β 1-42 ELISA Kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Target	A β 1-42
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cerebrospinal fluid (CSF)
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.