

## Product Datasheet

### Human IFN- $\gamma$ ELISA Kit (orb405942)

|                          |                                                        |
|--------------------------|--------------------------------------------------------|
| <b>Description</b>       | Human interferon gamma;(IFN-gamma;) antibody ELISA kit |
| <b>Reactivity</b>        | Human                                                  |
| <b>Range</b>             | Request Information                                    |
| <b>Target</b>            | IFN- $\gamma$ Ab                                       |
| <b>Note</b>              | For research use only                                  |
| <b>Application notes</b> | Detection Wavelength: 450 nm                           |
| <b>Sample Types</b>      | serum, plasma                                          |
| <b>Sensitivity</b>       | Request Information                                    |
| <b>Expiration Date</b>   | Please enquire.                                        |