

## Product Datasheet

### Human GGT2 ELISA Kit (orb405947)

<b>Description</b>	Human Gamma-glutamyltranspeptidase 2(GGT2) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	Request Information
<b>Target</b>	GGT2
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Uniprot ID</b>	<a href="#">P36268</a>
<b>Sensitivity</b>	Request Information
<b>Expiration Date</b>	Please enquire.