



Product Datasheet

RPL17 antibody (orb40595)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	RPL17
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	21397
Uniprot ID	P18621
NCBI	NP_001186274.1, NP_001186270.1, NP_001030178.1, NP_000976.1, NP_001186272.1, NP_001186271.1, NP_001186284.1, NP_001186269.1, NP_001186273.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of A375, HL-60 cel...