

## Product Datasheet

### Human Orexin A ELISA Kit (orb405951)

|                          |                                   |
|--------------------------|-----------------------------------|
| <b>Description</b>       | Human Orexin A ELISA Kit          |
| <b>Reactivity</b>        | Human                             |
| <b>Range</b>             | 0.156 ng/mL-10 ng/mL              |
| <b>Target</b>            | Orexin A                          |
| <b>Note</b>              | For research use only             |
| <b>Application notes</b> | Detection Wavelength: 450 nm      |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates |
| <b>Uniprot ID</b>        | <b>O43612</b>                     |
| <b>Sensitivity</b>       | 0.039 ng/mL                       |
| <b>Expiration Date</b>   | Please enquire.                   |