

Product Datasheet

Human CEACAM6 ELISA Kit (orb405954)

Description	Human Carcinoembryonic antigen-related cell adhesion molecule 6(CEACAM6) ELISA kit
Reactivity	Human
Range	0.83 ng/mL-20 ng/mL
Target	CEACAM6
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P40199
Sensitivity	0.33 ng/mL
Expiration Date	Please enquire.