

Product Datasheet

Human IL1 alpha ELISA Kit (orb405978)

Description	Human Interleukin 1 α , IL-1 α ELISA Kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	IL1A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
Uniprot ID	P01583
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.