

## Product Datasheet

### RGR antibody (orb40598)

**Description** Peptide affinity purified rabbit polyclonal antibody.

**Species/Host** Rabbit

**Reactivity** Human

**Conjugation** Unconjugated

**Tested Applications** FC, IHC-P, WB

**Immunogen** Synthetic Peptide

**Target** RGR

**Form/Appearance** Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

**Note** For research use only

**Isotype** Rabbit IgG

**Clonality** Polyclonal



**MW** 31874

**Uniprot ID** [P47804](#)


**NCBI** [NP\\_001012740.1](#), [NP\\_001012738.1](#), [NP\\_002912.2](#)

**Dilution Range** WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50

**Expiration Date** 12 months from date of receipt.



Immunohistochemical staining of paraffin...



Western blot analysis of in HL-60 cell l...