

## Product Datasheet

### Human IL1 Receptor ELISA Kit (orb405981)

<b>Description</b>	Human soluble interleukin-1 receptor I,IL-1sRI ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156 ng/ml-10 ng/ml
<b>Target</b>	IL1R1
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma,tissue homogenates
<b>Uniprot ID</b>	<a href="#">P14778</a>
<b>Sensitivity</b>	0.039 ng/ml
<b>Expiration Date</b>	Please enquire.