

Product Datasheet

Human IRAK4 ELISA Kit (orb405985)

Description	Human interleukin-1 receptor-associated kinase 4(IRAK-4)ELISA Kit
Reactivity	Human
Range	6.25 pg/mL-400 pg/mL
Target	IRAK4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9NWZ3
Sensitivity	1.56 pg/mL
Expiration Date	Please enquire.