

Product Datasheet

Human IL23A ELISA Kit (orb405994)

Description	Human Interleukin-23 subunit alpha(IL-23A) Elisa kit
Reactivity	Human
Range	62.5 pg/mL-4000 pg/mL
Target	IL23A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
Uniprot ID	Q9NPF7
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.