

Product Datasheet

Human IL26 ELISA Kit (orb405998)

Description	Human Interleukin 26,IL-26 ELISA Kit
Reactivity	Human
Range	62.5 pg/mL-4000 pg/mL
Target	IL26
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, urine, tissue homogenates, cell lysates, cell culture supernates
Uniprot ID	Q9NPH9
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.