

Product Datasheet

Human IL1F6 ELISA Kit (orb406011)

Description	Human Interleukin-36 alpha(IL1F6) ELISA kit
Reactivity	Human
Range	15.6 pg/mL-1000 pg/mL
Target	IL36A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell culture supernates
Uniprot ID	Q9UHA7
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.