

Product Datasheet

Human LILRB4 ELISA Kit (orb406027)

Description	human leukocyte immunoglobulin-like receptor subfamily B member 4, LILRB4 ELISA Kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	LILRB4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, cell lysates
Uniprot ID	Q8NHJ6
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.