

Product Datasheet

Human EFEMP1 ELISA Kit (orb406038)

Description	Human EGF-containing fibulin-like extracellular matrix protein 1(EFEMP1) ELISA kit
Reactivity	Human
Range	1.56 ng/mL-100 ng/mL
Target	EFEMP1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q12805
Sensitivity	0.39 ng/mL
Expiration Date	Please enquire.