

Product Datasheet

Human SOD3 ELISA Kit (orb406040)

Description	Human Extracellular superoxide dismutase [Cu-Zn](SOD3) ELISA kit
Reactivity	Human
Range	7.8 pg/mL-500 pg/mL
Target	SOD3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, saliva, urine
Uniprot ID	P08294
Sensitivity	1.95 pg/mL
Expiration Date	Please enquire.