

Product Datasheet

Human SEMA3A ELISA Kit (orb406044)

Description	Human Semaphorin 3A(Sema 3A)ELISA Kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	SEMA3A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q14563
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.