

Product Datasheet

Human HER2 ELISA Kit (orb406050)

Description	Human epidermal growth factor receptor 2 (sp185/HER2) ELISA Kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	ERBB2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P04626
Sensitivity	0.04 ng/mL
Expiration Date	Please enquire.