

Product Datasheet

Human EGFL7 ELISA Kit (orb406051)

Description Human Epidermal growth factor Like Domain Protein, Multiple 7 (EGFL7) ELISA Kit

Reactivity Human

Range 0.47 µg/mL-30 µg/mL

Target EGFL7

Note For research use only

Application notes Detection Wavelength: 450 nm

Sample Types serum, plasma, tissue homogenates

Uniprot ID **Q9UHF1**

Sensitivity 0.12 µg/mL

Expiration Date Please enquire.