

Product Datasheet

Human C1QTNF1 ELISA Kit (orb406072)

Description	Human Complement C1q tumor necrosis factor-related protein 1(C1QTNF1/CTRP1) ELISA Kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Target	C1QTNF1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9BXJ1
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.