

Product Datasheet

Human C1QTNF2 ELISA Kit (orb406073)

Description Human Complement C1q tumor necrosis factor-related protein 2(C1QTNF2) ELISA kit

Reactivity Human

Range 0.156 ng/mL-10 ng/mL

Target C1QTNF2

Note For research use only

Application notes Detection Wavelength: 450 nm

Sample Types serum, plasma, tissue homogenates

Uniprot ID **Q9BXJ5**

Sensitivity 0.039 ng/mL

Expiration Date Please enquire.