

## Product Datasheet

### Human MCP1 ELISA Kit (orb406111)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | Human high sensitivity C-C motif chemokine 2 (CCL2) ELISA kit |
| <b>Reactivity</b>        | Human   |
| <b>Range</b>             | 6.25 pg/mL-400 pg/mL  |
| <b>Target</b>            | CCL2  |
| <b>Note</b>              | For research use only   |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                  |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                             |
| <b>Uniprot ID</b>        | <b>P13500</b>   |
| <b>Sensitivity</b>       | 1.56 pg/mL  |
| <b>Expiration Date</b>   | Please enquire.   |