

Product Datasheet

Human hs-CRP ELISA Kit (orb406113)

Description	Human high sensitivity C-Reactive Protein,hs-CRP ELISA Kit
Reactivity	Human
Range	0.625 ng/mL-40 ng/mL
Target	hs-CRP
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P02741
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.