

Product Datasheet

Human FGF23 ELISA Kit (orb406118)

Description	Human fibroblast growth factor-23,FGF-23 ELISA Kit
Reactivity	Human
Range	3.12 pg/mL-200 pg/mL
Target	FGF23
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, cell culture supernates, urine, tissue homogenates
Uniprot ID	Q9GZV9
Sensitivity	0.78 pg/mL
Expiration Date	Please enquire.