



Product Datasheet

ICB1 antibody (orb40612)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	THEMIS2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	72049
Uniprot ID	Q5TEJ8
NCBI	NP_001034566.1 , NP_001099026.1 , NP_004839.2
Dilution Range	WB 1:50-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of in HL-60 cell l...